Smart steps for first-time wearers

It takes practice — and patience to get comfortable with hearing aids. These quick tips should help speed up the process:

- Wear them at home in a quiet environment for the first few days
- Practice having a conversation with your spouse or a loved one
- Read out loud for 10-15 minutes a day
- Wear them only a few hours each day for the first two weeks
- **** Take breaks when you feel tired
- Do online hearing skills training exercises (your hearing professional can provide fun interactive games)
- Set realistic expectations; hearing aids won't restore your hearing to normal — but they will make listening easier and much more enjoyable again




Starkey Hearing Technologies 6700 Washington Avenue South Eden Prairie, MN 55344

© 2014 Starkey Hearing Technologies. All Rights Reserved. 85404-007 10/14 BR0C0324-02-EE-SG

What to Expect From Your Hearing Aids

How they help, and what you can do to get the most out of them.





change your perception and get ready to change your life

Most people don't know what to expect from hearing aids. Misperceptions and second-hand experience with bulky, whistling, old-fashioned analog devices continue to influence the way people think about all hearing aids.

Hearing aids have changed dramatically from a generation ago — from outward appearance to internal technology making them vastly more appealing and effective.

This brochure will help you set expectations for the performance of your new hearing aids — how they sound, feel and work.

how your hearing aids Should Getting acclimated to wearing your new

Getting acclimated to wearing your new hearing aids is different for each person (sometimes taking up to two to three months). But once you are, you should notice a big difference in how the world sounds.

- Your ability to hear and understand others should be improved
- Other people's voices shouldn't sound distorted, harsh, tinny, sharp, booming or muffled
- The sound of your own voice should be "normal," not sound like you're in a barrel
- The intensity and quality of familiar sounds should be sharp, bright and clear not dull or irritating
- In a crowded room with many people talking at once, sound should be loud but not deafening
- Wearing hearing aids in both ears should help you identify the location of a sound or voice
- Hearing and communicating in quiet environments (home, work, doctor's office) should be improved
- Your ability to hear and understand speech in environments with background noise (restaurants or dinner parties, for instance) should be improved
- Your hearing aids should help you understand speech in larger environments where there is reverberation (lecture halls, worship spaces, movie theaters)
- Loud sounds (sirens, traffic, construction sounds) should not be uncomfortable but you should hear them clearly

Five tips to get even more out of your hearing aids:

- 1 If you wear directional hearing aids, sit facing the wall with the restaurant noise behind you
- 2 Take them into your hearing professional for regular cleaning and maintenance
- 3 Make sure to keep fresh batteries on hand
- 4 Turn off hearing aids or remove battery when not in use
- 5 In public places with acoustic challenges (theaters or conference halls), it is best to sit in the front and center of the room, where it offers the best acoustics

If you have any questions about your hearing aids, talk to your hearing healthcare professional.