SONIC captivate.



miniRITE TR in Charger

Some Captivate devices are rechargeable, so you can enjoy all-day power, even when streaming. Ask your Hearing Care Provider about models and details.

Bluetooth is a registered trademark owned by Bluetooth SIG, Inc, USA.

Apple, the Apple logo, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Android and Google Play are trademarks of Google LLC.

© 2019 Sonic Innovations, Inc. All rights reserved.

2019-01-24 | 203952 | UK



What life sounds like

Enjoy more of the world of sound, with less noise and fewer distractions. From clarifying conversations to streaming online entertainment, Captivate, an advanced new hearing solution, is here today to immerse you in the sounds of life.



Captivate miniRITE TR

Built-in rechargeability

≰iPhone | iPad | iPod

Experience the sounds of life

Introducing Captivate—better hearing throughout the moments of your life.*

Captivate uses the latest innovations in the Sonic **SoundDNA** technology platform to keep conversations clear, even in noisy environments. Feedback is virtually eliminated, keeping distractions to a minimum. You can easily connect to the online world, to your smartphone and other wireless devices.

Small, stylish and highly customizable, Captivate is resistant to dust and moisture, and can even provide a full day of power with optional rechargeability.

Captivate makes everyday sounds better in exciting new ways.*

Natural sound

SmartCompress analyzes your listening environment and adjusts to sound quickly and automatically for a natural sound experience.

Clear speech

Speech in Noise (SPIN) Management improves your ability to hear words in noisy environments, where you may need the most help.

Fewer distractions

NEW! Adaptive Feedback Canceller Pro quickly removes feedback and helps keep your ear canal open, so you can hear the sound of your own voice more naturally compared to our previous technology.

Inspiring experiences

Extended Dynamic Range, available only in Captivate¹⁰⁰, keeps loud, dramatic sounds clear, not distorted, in places like movie theaters and auditoriums, for more impactful entertainment.



SoundDNA Platform

Personalize Captivate for a unique hearing experience.*

Power through the day

NEW! Lithium-ion rechargeability lasts all day, even while streaming. With fast, efficient charging, this built-in feature on the **miniRITE T R** is available in three technology levels.

Ease into your new hearing aids

Captivate's **Adaptation Manager** automatically adjusts your settings over a period of time until you're comfortable with the volume and clarity you need.

Reduce the ringing in your ears

Built-in **Tinnitus SoundSupport** generates relief sounds suited to your preferences to reduce the perception of tinnitus.

Help for high-frequency loss

If you have a severe high-frequency hearing loss, certain high-pitched sounds, like women's and children's voices, may be more difficult to hear. **Frequency Transfer** shifts these sounds to a range that's easier for you to perceive.

* Not all features and options available on all technology levels or models. Ask your Hearing Care Provider for details.

Wonderfully Wireless



Our **Dual-Radio System** uses 2.4 GHz Bluetooth® Low Energy (BLE), so **no intermediary device is needed** to stream sound from select digital sources to your hearing aids. Also included is Near Field Magnetic Induction (NFMI) technology, so information is traded between right and left hearing aids quickly for a natural sound experience.



Captivate is a Made for iPhone® hearing aid, compatible with devices running iOS 9.3 or later. For information on compatibility, please visit www.sonici.global/compatibility.

A variety of optional wireless accessories can help you create new listening experiences, every day.

SoundLink 2 App-Change your listening program, adjust volume, control tinnitus relief sounds and switch between accessories conveniently from your iPhone® or Android™ smartphone.¹





SoundClip-A–Enjoy stereo streaming with Bluetooth-compatible devices. Also functions as a remote microphone, remote control and wireless headphones for hands-free calls. Use **Phone Adapter 2** with SoundClip-A to stream calls from your landline phone directly to Captivate.

RC-A Remote Control-Adjust volume and programs with a click. Small enough for pocket or purse.

TV-A Adapter–Stream audio from your TV directly to Captivate hearing aids. This adapter connects easily to most TVs.

If This Then That (IFTTT)—Connect your hearing aids to the Internet of Things with IFTTT, a web-based service that can integrate sound cues into your life. For example, you can receive an audible signal when you get a new email or a different sound alert when someone's at the front door. Become more aware of the world around you. Requires SoundLink 2 App.

Express your style

Captivate is now available in three technology levels and four elegant models, with up to seven sophisticated colors to choose from. Ask your Hearing Care Provider which combination works best for you.

For more information about Captivate, visit www.sonici.com.



* Gold color available in the Captivate 100 miniRITE only.