


ReSound ONE

Made for **É** iPhone | iPad | iPod

Hear like no other



Introducing ReSound ONE

Just like a fingerprint, each person's hearing is ONE of a kind. This is because nature has designed the shape of your ear to perfectly capture the sounds around you.

ReSound ONE™ with M&RIE (Microphone and Receiver-in-Ear) delivers a new class* of hearing aid with an additional microphone inside your ear, allowing the unique shape of your ear to collect sounds, just as nature intended.

Add convenient direct audio streaming, access to remote hearing care services, and industry-leading battery life, ReSound ONE truly keeps you connected to people and technology you need to become ONE with your world.



 ${\rm *Microphone}\ in\ ear\ and\ two\ standard\ directional\ microphones, directionality\ features, wireless\ audio\ streaming.$



ReSound ONE provides you with a more complete and individualized sound experience with:

A MORE NATURAL SOUND AND IMPROVED SENSE OF SPACE

ReSound ONE with M&RIE is the first hearing aid* with three microphones, including one in the ear, to use your ear's natural ability to collect the information your brain needs to make sense of sound.

EASIER-TO-FOLLOW CONVERSATIONS IN DYNAMIC ENVIRONMENTS

All Access Directionality gives you access to all the sounds around you, putting you in charge of which sounds you want to focus on - not your hearing aids.

BEST 1:1 SPEECH UNDERSTANDING

Like an invisible beam of sound detection, Ultra Focus lets you decide when to concentrate on the speech in front of you while minimizing background noise - even in the toughest environments.

Unmatched connectivity to adapt to your needs

ReSound ONE adapts to your lifestyle connecting you to the world around you with direct audio streaming from both iOS^{\otimes} and AndroidTM devices. So your favorite audio book, music or even GPS directions can stream directly to your hearing aids without an intermediary device.

One-tap sound adjustments on your ReSound Smart 3D™ app and a versatile range of wireless accessories give you even more options to connect to the sounds you enjoy every day, like your favorite TV program.

More connections for better care


ReSound ONE includes access to ReSound Assist, so you can receive hearing care in the most convenient way.

FACE-TO-FACE ONLINE HEARING CARE

ReSound Assist Live Assistance makes it possible to have face-to-face fitting sessions from the comfort of your home. Simply answer the incoming video call from your hearing care professional and get real-time fine-tuning adjustments.

REMOTE FINE-TUNING ANYTIME

ReSound Assist Remote Fine-Tuning allows you to send requests for additional adjustments anytime through the ReSound Smart 3D app. Updates will be sent to you to download at your convenience - without the need for an office visit.



Designed for ultimate comfort, discretion and performance

Built for ultimate comfort and adaptability, the modern styling of ReSound ONE has just one touchpoint on the back of the ear and sits lower for greater discretion. Available in three models, ReSound ONE comes in a variety of colors to blend with hair and skin tones, so you can wear ReSound ONE with confidence.



INDUSTRY-LEADING RECHARGEABILITY

ReSound ONE offers the best rechargeable battery life in the industry with up to 30 hours* of battery life or 25 hours with the ability to stream as much as you want.

For added flexibility, two different charger options serve as stylish, portable and protective cases, ideal for traveling and storage. Choose the Premium Charger with an onboard battery for up to three days of quick charging on the go, without the need for a power outlet, or the smaller Standard Charger that gives you at at-home charging capabilities you need in a small elegant design.





At ReSound we believe life is best when it's a natural, effortless flow, so we take to heart anything that holds you back. This is why we have devoted the past 150 years to making life sound better, connecting you to a life without limitations. Our breakthrough innovations, inspired by our Organic Hearing philosophy, all have one core ingredient - you, because hearing and hearing care should be as individual as a fingerprint.

ReSound is part of the GN Group - pioneering great sound from world-leading ReSound hearing aids to Jabra office headsets and sports headphones. The GN Group was founded in 1869, employs over 5,000 people. and is listed on NASDAQ OMX Copenhagen.

Find out how you can hear like no other with ReSound ONE.

resound.com

facebook.com/resoundhearing

twitter.com/resoundus

youtube.com/resoundus

Manufacturer according to FDA:

GN ReSound North America

8001 E Bloomington Freeway Bloomington, MN 55420 USA 1-888-735-4327

resound.com

Manufacturer according to Health Canada:

ReSound Canada

2 East Beaver Creek Road, Building 3 Richmond Hill, ON L4B 2N3 Canada 1-888-737-6863 resound.com

© 2020 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google Inc.